



WEIBER[®]



SOXHLET EXTRACTION UNIT

WEIBER[®]

Models:

- ACM-54096 W



ACMAS TECHNOCRACY PVT. LTD.

www.measuring-meters.com | www.acmasindia.com | www.scientificalaboratoryequipments.com

WEIBER[®]

Process and Technology Explanation ;

The equipment was originally designed for the extraction of a lipid from a solid material. However, a Soxhlet extractor is not limited to the extraction of lipids. Typically, a Soxhlet extraction is only required where the desired compound has a limited solubility in a solvent, and the impurity is insoluble in that solvent. If the desired compound has a high solubility in a solvent then a simple filtration can be used to separate the compound from the insoluble substance.

Normally a solid material containing some of the desired compound is placed inside a thimble made from thick filter paper, which is loaded into the main chamber of the Soxhlet extractor. The Soxhlet extractor is placed onto a flask containing the extraction solvent. The Soxhlet is then equipped with a condenser.

The solvent is heated to reflux. The solvent vapour travels up a distillation arm, and floods into the chamber housing the thimble of solid. The condenser ensures that any solvent vapour cools, and drips back down into the chamber housing the solid material. The chamber containing the solid material slowly fills with warm solvent. Some of the desired compound will then dissolve in the warm solvent. When the Soxhlet chamber is

SOXHLET EXTRACTION UNIT

ACMAS TECHNOCRACY PVT. LTD.

www.measuring-meters.com | www.acmasindia.com | www.scientificalaboratoryequipments.com

The logo for WEIBER, featuring the brand name in a bold, black, sans-serif font with a registered trademark symbol (®) to the upper right of the 'R'. The logo is centered between two horizontal grey bars.The logo for WEIBER, featuring the brand name in a bold, black, sans-serif font with a registered trademark symbol (®) to the upper right of the 'R'. The logo is centered between two horizontal grey bars.

SOXHLET EXTRACTION UNIT

almost full, the chamber is automatically emptied by a siphon side arm, with the solvent running back down to the distillation flask. This cycle may be allowed to repeat many times, over hours or days.

During each cycle, a portion of the non-volatile compound dissolves in the solvent. After many cycles the desired compound is concentrated in the distillation flask. The advantage of this system is that instead of many portions of warm solvent being passed through the sample, just one batch of solvent is recycled.

After extraction the solvent is removed, typically by means of a rotary evaporator, yielding the extracted compound. **The non-soluble** portion of the extracted solid remains in the thimble, and is usually discarded.

Product Construction Details

The complete soxhlet extraction apparatus consists of solvent reservoir extraction assembly with 6 extraction places with concentric metal rings of various diameter to suit 100 ml flasks to

SOXHLET EXTRACTION UNIT

250 ml flasks and built in steam generator. Complete with stand, clamp holder soxhlet extraction (Condenser) glass and extraction flask of 60 ml volume.

Temperature Range :

Temperature Range : 50°C to 350°C

Temperature Control:

Temperature is controlled through electronic solid state temperature controller cum indicator or micro processor based digital temperature controller cum indicator or hydraulic thermostat with energy regulators.

Temperature Sensitivity:

The temperature is controlled with an efficiency of $\pm 2\%$ (Set Value)

Power Requirements:

Power 2 KW, Single Phase, 230 Volts.



We are having **Representatives in 55 Countries**


please visit our

International Dealers Section

to contact your nearest **ACMAS** Representative



TECHNOLOGY with **HUMAN TOUCH**



CE | An ISO 9001:2208 | ISO 14001:2008 | ISO 13485
WHO:GMP Products | GLP Compliant Products

www.acmasindia.com | www.scientificlaboratoryequipments.com | www.measuring-meters.com

ACMAS TECHNOCRACY PVT. LTD.

SALES OFFICE (INDIA)

312-313, Vardhman Capital Mall, L.S.C. 10, GulabiBagh, Delhi-110052, INDIA

Tel: 0091-11-23646703, 23643054

(M) +91-9313971681 | Email: india@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/f Wing Tat Comm, Bldg 97, Bonham Strand East, Sheung Wan, HONG KONG (PRC).

Tel: 0086-13929598046 | 0086-18922303099

Email: hk@acmasindia.com

SHOWROOM

141, Rai Industrial Estate, Rai, Sonapat, Haryana-131029, INDIA

Tel: 0091-0-9312219738 (M) +91-9717741167

Email: info@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street, Ozerkovsky Embankment, Unit No 50 Straine-1, Off-502, Moscow, RUSSIA

Tel: 0049-79592345 | Email: russia@acmasindia.com