

# REFRIGERATED CENTRIFUGE

Multipurpose Centrifuge that accommodates a large number of Rotor Heads of varying tube capacities. The Table-top model centrifuge is provided with step-less autovariac speed control and with separate control panel having digital speed indicator, digital Timer of 0-99 minutes and automatic switching off, alarm signal, digital temperature control cum indicator with a range of -20 oC to +40 o c. The equipment also feature dynamic barke, zero start inter lock, safety cut off in case of imbalance and lid locking switch.



## Specifications

<b>Acm 67894-C</b>			
Max. Speed	RPM		18000
Max. RCF	x g		28600
Max. Tube Size	ml		200
Max. Capacity	ml		800
Lowest Temp.	°C		-10
Width	mm		875
Depth	mm		620
Height	mm		660
Weight	kg		120
Connected Load	KVA		1.75

## Interchangeable Rotor Head Suitable For Use With ACM 67894-C Centrifuge

<b>Model</b>	<b>Cap. / ml</b>	<b>RPM</b>	<b>RCF x g</b>
AC 8101 A	4 x 200	5200	5050
AC 8101 B	4 x 100	6100	6600
AC 8101 C	4 x 18 RIA	5200	4200
AC 8101 D + E	4 x 27 RIA	5200	4200
<b>OTHER AVAILABLE SWING OUT ROTOR HEADS</b>			
AC 8107	2 x Microtiter Trays	3400	1700
AC 4110	4 x 10	11000	12900
AC 4114	4 x 25	11000	13800

# REFRIGERATED CENTRIFUGE

## FIXED ANGLE ROTOR HEAD

<b>Model</b>	<b>Cap. / ml</b>	<b>RPM</b>	<b>RCF x g</b>
AC 8102	12 x 25	7000	7100
AC 8103	6 x 200	7000	7250
AC 8104	8 x 100	7000	7100
AC 8105	12 x 50	7000	7100
AC 8106	36 x 15	7000	7100
AC 4101	4 x 100	14000	21800
AC 4102	8 x 50	13500	21800
AC 4103	6 x 50	14500	23000
TC 4104	4 x 50	16000	23500
TC 4105	8 x 25	15000	23000
TC 4106	6 x 25	16000	25200
TC 4107	12 x 15	15000	23000
TC 4108	8 x 15	17500	28600
TC 4109	6 x 15	18000	27900
TC 4111	24 x 1.5	18000	25500
TC 4112	24 x 1	18000	25500
TC 4113	24 x 0.5	18000	25500

# PROFILE

## COMPANY PROFILE :

With over a decade of industry experience in manufacturing scientific laboratory instruments, we have established ourselves with the name of **ACMAS Technocracy (P) Ltd.** to provide highest quality instruments to laboratories, pathologies, entomologies, pharmaceuticals, research centers etc. We have successfully catered the needs of above 600 institutions in India and abroad for the last 22 years. The dedicated and cumulative efforts of ACMAS team members has produced and delivered the comprehensive range of scientific instruments and laboratory products research projects where maintaining a viable record of the performance of the equipment is very essential.

## COMPANY'S OUTLOOK :

**ACMAS Technocracy (P) Ltd.** enjoys an amazing image for high quality scientific laboratory instruments across the globe. The continuous innovative technology and 'Quality Management System Standard' delivers the advanced laboratory experiments and general-purpose measuring instruments solution to various laboratories, sterilizing clean rooms, microbiologies, biotechnologies, pathologies, entomologies, pharmaceuticals, seeds and soil testing, meteorologies food processing. We also believe in providing customized instruments solution to our esteemed clients.

## MANUFACTURING FACILITIES :

The company has well built and operated manufacturing facilities that matches the latest system and technique in the industry. Our team strictly follows the quality control standards of ISO 9001:2000 series while designing, developing, manufacturing and delivering the scientific instruments. The manufacturing unit of **ACMAS** is made with complete state-of-the-art equipments and technologies for producing high quality instruments. We also acquired Environmental Friendly process certifications ISO14001:2004 for our entire range of instruments to ensure reliability and durability in each product.

## RANGE OF PRODUCTS :

Our wide array of **AUTOClave, INCUBATORS OVENS, LAMINAR AIR FLOW, MOISTURE METERS, WATER DISTILLATION PLANTS, LABORATORY BALANCES, WATER BATH, CENTRIFUGE, COOLING EQUIPMENTS, WATER TESTING EQUIPMENTS, LABORATORY SHAKING MACHINE, MICROTOME, MICROSCOPE, MEASURING INSTRUMENTS AND ALLIED** products to ensures accuracy and conformity for significant experiments. We also customize some of our product range and technologies for educational, medical, industrial or other laboratories for better work experience. Our pre sale and post sales support are also admirable and popular among our satisfied clients.

## QUALITY STANDARDS :

At **ACMAS Technocracy (P) Ltd.**, we design and develop the complete range of scientific and laboratory instruments with highest quality standards. We constantly update technologies and methodologies to ensure reliability and consistency at each level of instruments production. Our all transparency auditing system are performed by the most reliable D& B International as we want to deliver the world class quality instruments to our national and international clients. We feel proud that our entire product range has brought satisfactory results for the corporate and public sector clients.

## PROFESSIONAL TEAM :

The continuous cooperation and support of professional team has helped us to understand and deliver the satisfactory scientific instruments right from basic lab equipment to most sophisticated instruments for research labs. We believe that our tremendous success belongs to our expert engineers, managers, co-workers and other significant team members who have put their best efforts in the growth of the organization. It is their dedication and commitment that makes us the most trusted scientific instruments brand among our all satisfied clients.

## FORTE :

We strive hard to cater our clients with best product range and services while meeting international standards. Our aim is to meet the overwhelming demand of the scientific community and provide them the world class quality scientific instruments along with the best after sale support.



# CLIENTS

OUR LIST OF CUSTOMERS :



# CLIENTS

OUR LIST OF CUSTOMERS :



# OUR PRESENCE WORLD WIDE



# CONTACT

## ACMAS TECHNOCRACY (P) LTD

### Marketing & Sales :

Atul Badola (Manager Marketing)

atulbadola@gmail.com, atul@acmasindia.com

### Production :

K.K. Pawar (Manager Production)

kk@acmasindia.com

### Technical Support :

Sandeep Bose (Manager Technical Support)

sandeep@acmasindia.com

### Pre and Post Sales Customer Support :

Sonia Nathani (Co-ordinator Customer Support)

sonia@acmasindia.com

### Business Development :

Abha Verma (Executive- Business Development)

abha@acmasindia.com

### Accounts and Finance :

Meghna Arora (Executive – Accounts)

meghna@acmasindia.com

### Product Development :

Nishu Tomar (Executive – Product Development)

nishu@acmasindia.com

## PHONE :

### Hand Phones :

+91-0-9717741167, +91-0-9312219738,  
+91-0-9313971681, +91-0-9350565689

### Land Line - Office :

+91-011-23643054, +91-011-23646703

### Land Line - Works :

+91-011-22942133, +91-011-22943508

### Telefax :

+91-011-23646703

## WORKS :

### Unit I :

A-100/1 Main Som Bazar Road, Gamri  
Extension PO: Maujpur,

Delhi - 110052 (INDIA)

### Unit II :

1/6 DSIDC Complex, Nand Nagri,

Delhi - 110 032 (INDIA)

## SALES OFFICE :

# 312-313 Vardhman Capital Mall,

Local shopping Complex, Gulabi Bagh

Delhi - 110054 (INDIA)

