






Models:  
ACM-78095 S



**WEIBER**

**Ice Flaker Machine**



Models:



- ACM-42103 L



# WEIBER




Ice Flaker Machine



# WEIBER

Ice Flaker Machine



**Indicative photographs only**

We are manufacturers of ice flaker machines in India since last two decades.. Our ice flaker machines are conceptualized and designed to meet the challenging to all the microbiology, biochemistry , biotechnology and genetic engineering laboratories that have a continuous requirement of ice to store their biological samples along with transportation to other facilities.

We strongly believe in the principle of quality first and everything else to follow, which reflects in our products, which are renowned for long trouble free services in various research and educational institutes across India and abroad.



# WEIBER



Ice Flaker Machine

Our ice flaker machines are made of corrosion free metal components with advance poly uretherane insulation ensuring maximum thermal efficiency and durability.

We specialize in both standard and customized models of ice flaker machines specially designed to meet various special requirements of our customers across India and abroad..

Over a short period of time Weiber brand have been established as a leading exporters of ice flaker machines in India, catering to South Asian, African And European Markets. Beside the regular exports we are also catering to our domestic market and are suppliers of our ice flaker machines to various blood banks, hospitals, research stations, R & D Laboratories educational institutes and industries.

#### Salient Features:

- Corrosion free exterior and interior
- Designed for continous production of ice throughout the day
- Environmental friendly CFC free cooling
- Micro processor based temperature control
- PUF insulation.
- Reliable
- Energy Efficient
- Sturdy Construction
- Low Maintenance
- Light Weight

ACMAS TECHNOCRACY PVT. LTD.

[www.measuring-meters.com](http://www.measuring-meters.com) | [www.acmasindia.com](http://www.acmasindia.com) | [www.scientificlaboratoryequipments.com](http://www.scientificlaboratoryequipments.com)



# WEIBER



Ice Flaker Machine

#### Technical Features:

- Maximum Production : 50-100 kg / Day (During 24 hour operation)
- Quality compressor & Air Cooled condensation CFC/HCFC Free. (R134a)
- Storage Capacity of bin : 16 kg
- Type of ice : flakes/Grains
- Interior and Exterior are made of corrosion resistant material.
- Easy Extraction system
- Energy Efficient
- Automatic low water cut off system
- Automatic full storage cut off system
- Provided with geared motor with over heating and over voltage protection
- Production start time : 10-15 mins
- Provided with castor wheels for easy movement

ACMAS TECHNOCRACY PVT. LTD.

[www.measuring-meters.com](http://www.measuring-meters.com) | [www.acmasindia.com](http://www.acmasindia.com) | [www.scientificlaboratoryequipments.com](http://www.scientificlaboratoryequipments.com)



We are having **Representatives in 55 Countries**

please visit our

**International Dealers Section**

to contact your nearest **ACMAS** Representative



**TECHNOLOGY** with **HUMAN TOUCH**



An ISO 9001:2208 | ISO 14001:2008 | ISO 13485  
WHO:GMP Products | GLP Compliant Products

[www.acmasindia.com](http://www.acmasindia.com) | [www.scientificlaboratoryequipments.com](http://www.scientificlaboratoryequipments.com) | [www.measuring-meters.com](http://www.measuring-meters.com)

## **ACMAS** TECHNOCRACY PVT. LTD.

### **SALES OFFICE (INDIA)**

312-313, Vardhman Capital Mall, L.S.C. 10, GulabiBagh, Delhi-110052, INDIA

Tel: 0091-11-23646703, 23643054

(M) +91-9313971681 | Email: [india@acmasindia.com](mailto:india@acmasindia.com)

### **SALES OFFICE (HONG KONG)**

Unit D 28 11/f Wing Tat Comm, Bldg 97, Bonham Strand East, Sheung Wan, HONG KONG (PRC).

Tel: 0086-13929598046 | 0086-18922303099

Email: [hk@acmasindia.com](mailto:hk@acmasindia.com)

### **SHOWROOM**

141, Rai Industrial Estate, Rai, Sonapat, Haryana-131029, INDIA

Tel: 0091-0-9312219738 (M) +91-9717741167

Email: [info@acmasindia.com](mailto:info@acmasindia.com)

### **SALES OFFICE (RUSSIA)**

Inmed Trade Street, Ozerkovsky Embankment, Unit No 50 Straine-1, Off-502, Moscow, RUSSIA

Tel: 0049-79592345 | Email: [russia@acmasindia.com](mailto:russia@acmasindia.com)