



Weiber Ultra Low Deep Freezer

MODEL NO. - (ACM- DW-HL) 218/328/388/538/828)





Weiber Ultra Low Deep Freezer

MODEL NO. - (ACM- DW-HL) 218/328/388/538/828)

PRODUCT DESCRIPTION

Effective Capacity: 218/328/388/538/828 Liter

Temperature Range: 10°C~ -86°C (limit temp. is -90°C)

FEATURES

- Upright type
- Inside temperature ranging from -10°C~-86°C can be set freely , controlling precision is 1°C
- Stainless steel interior, Outer is painted steel panel
- Front opening lockable door with full height handle
- 155mm foaming insulation (DW-HL538 is 100mm with VIP board)
- Vacuum release port
- Inner 3 shelves Stainless steel with inner door (DW-HL218 is 2shelves)
- Imported Danfoss compressor, Tow compressors parallel connection, to reduce noise and improve the cooling temperature come to -86°C steady.
- 2-4 hours come to -86°C .very fast.
- Fan electromotor with internationally famous of Germany EBM



- Microprocessor-based temperature control system, platinum resistance sensors
- Digital temperature display
- Keyboard lock and password protected configuration page
- Audible and visual alarm: High or low temperature alarm, Power failure alarm, Low battery alarm, Door ajar alarm, Filter blocking alarm, system failure alarm
- With temperature printer, record 7days data
- 72 hour battery backup for printer and controller display
- 4units Casters for easy handing
- Mixture gas with R134a,R404a,R23,R24,N2,R600a and more, CFC free
- Certificate: CE mark for DW-HL328/538, ISO9001, ISO14001,ISO13485
- Power supply: 220V /50Hz 1 phase, can change as 220V 60HZ or 110V 50HZ or 110V 60HZ

OPTIONAL: CHART RECORDER, CO2 BACK UP, STORAGE RACKS, BOXES

Model	Capacity	External size (WXDXH)mm	Inside size (WxDxH) mm	Input power	Weight (N.t / G.t)
DW-HL218	218 liters	940x840x1552	582x470x768	1000W	225/250KG
DW-HL328 (CE mark)	328 liters	940x840x1999	582x470x1200	1500W	320/350KG
DW-HL388	388 liters	940x840x1999	600x538x1200	1800W	330/360KG
DW-HL538(CE mark)	538 liters	940x840x1999	650x600x1345	2000W	360/380KG
DW-HL828	828 liters	1230x996x2056	930x710x1323	3200W	420/450KG



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485
WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.acmasglobal.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57
Phase- IV Kundli, Sonapat, HR 131028
Land Line No.: +91 - 11- 47619688,
Fax: +91-11-47619788
E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,
Sheung Wan, Hong Kong (PRC)
Tel.: 0086-13929598046 | 0086- 18922303099
E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,
Unit No 50, Straine- 1, Off- 502, Moscow, Russia
E-mail.: info@acmasindia.com
Tel.: 0049- 79592345 | Email: russia@acmasindia.com