



DEEP FREEZER

Weiber Deep Freezer

MODEL No. : ACM-DF-4559



WEIBER DEEP FREEZER

MODEL No.: ACM-DF-4559



Vertical Deep Freezer



Horizontal Deep Freezer

FEATURES

- * Compact design with large storage capacity
- * Outer casing is galvanized steel with baked on paint finish for durability.
- * Robust construction and elegant looks.
- * Inner liner is Stainless Steel/Pre coated GI
- * Higher Capacity System for fast pulls down.
- * High density, environment friendly, polyurethane foam for better insulation.
- * High-pressure injection machine used for efficient insulation.
- * R-134 a refrigerant, chlorine free, environment friendly.
- * Carrier/Kirloskar/Tecumseh compressor.
- * Built in voltage stabilizer with time delay relay
- * Energy saver compressor-low energy consumption.
- * Longer thaw out time to take care of prolonged periods of failure.
- * 60mm thick insulation
- * Temperature range (ambient to -22°C)
- * Temperature controlled through solid-state digital temperature controller cum indicator.
- * Volume: 250 liters / 350 Liters / 500 Liters (It can be customized as per the requirement)

Optional: Voltage Stabilizer (2 KVA) Auto Cut Off Stabilizer

Note: We can also customized the product as per the customer requirements.





An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485
WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.acmasglobal.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57
Phase- IV Kundli, Sonapat, HR 131028
Land Line No.: +91 - 11- 47619688,
Fax: +91-11-47619788
E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,
Sheung Wan, Hong Kong (PRC)
Tel.: 0086-13929598046 | 0086- 18922303099
E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,
Unit No 50, Straine- 1, Off- 502, Moscow, Russia
E-mail.: info@acmasindia.com
Tel.: 0049- 79592345 | Email: russia@acmasindia.com