

WEIBER[®]

Models:
ACM-78093-L



Multi-functional Infrared Moisture Analyzer

Models:

- ACM-78093-L

WEIBER[®]



Product Information

Our Multi-functional Infrared Moisture Analyzer ((Acm-78093 L)) comprises of a high precision analytical balance and a structured stepless heating device. This equipment is based on the principle of electromagnetic balancing principle combined with the infra red heating technology used at tandem to ascertain the moisture content in the sample.

Multi-functional Infrared Moisture Analyzer

This unique infra red moisture balance is powered by the integrated advanced micro processor, thus making it highly accurate and reliable. Beside the accuracy durability sturdiness and diligence are other key attributes of our products. This equipment is very user friendly allowing the user to either manually or automatically set the time temperature, mode and intensity of heating and other parameters in the unit.

The unit is provided with RS-232 C interface with computer connectivity . The user can have four output modes of selection with eight parameter output.

These features makes it an ideal choice for tobacco, pharma , food processing and chemical laboratories.

Multi-functional Infrared Moisture Analyzer

Technical Specifications

Features

- Advanced Infra Red Heating (Halogen)
- Large Liquid Crystal display
- User Friendly
- High Accuracy
- Good Repeatability
- Computer Interface Through RS 232 C Interface.

Application:

- Research Laboratories.
- Product Development Applications
- Product Testing Applications
- Pharmaceutical Testing
- Various Industrial Applications

Model	(Acm-78093 L)
Accuracy (Working)	0.001g
Cap	0-100g
Operating Temp Range	Ambient to 160°C
Accuracy (Temp Control)	±1°C
Readability (Moisture)	0.01%
Repeatability(Sample)	0.2%(sample ³ 2g)
Setting of time	0-120 min
Mode	Manual/ Automatic
Control Module	Program control
Display	8 Parameter
Dia (Pan)	90mm
Dimensions	325mm x 195mm x 195mm
Power Supply	220V+22V (50Hz) 110V±11V (60Hz) 220V-33V (50HZ)
Net Weight	5.8 kg



Multi-functional Infrared Moisture Analyzer



We are having **Representatives in 55 Countries**

please visit our

International Dealers Section

to contact your nearest **ACMAS** Representative



TECHNOLOGY with **HUMAN TOUCH**



An ISO 9001:2208 | ISO 14001:2008 | ISO 13485
WHO:GMP Products | GLP Compliant Products

www.acmasindia.com | www.scientificlaboratoryequipments.com | www.measuring-meters.com

ACMAS TECHNOCRACY PVT. LTD.

SALES OFFICE (INDIA)

312-313, Vardhman Capital Mall, L.S.C. 10, GulabiBagh, Delhi-110052, INDIA

Tel: 0091-11-23646703, 23643054

(M) +91-9313971681 | Email: india@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/f Wing Tat Comm, Bldg 97, Bonham Strand East, Sheung Wan, HONG KONG (PRC).

Tel: 0086-13929598046 | 0086-18922303099

Email: hk@acmasindia.com

SHOWROOM

141, Rai Industrial Estate, Rai, Sonapat, Haryana-131029, INDIA

Tel: 0091-0-9312219738 (M) +91-9717741167

Email: info@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street, Ozerkovsky Embankment, Unit No 50 Straine-1, Off-502, Moscow, RUSSIA

Tel: 0049-79592345 | Email: russia@acmasindia.com