



GENERAL PURPOSE PHYSICAL **LABORATORY BALANCE**

MODEL NO. - ACM- 78091- L



General Purpose **Physical Laboratory** **Balance**

FEATURES

- Liquid Crystal Display
- High Internal Resolution
- Grams & Carats.
- Digital Auto Calibration.
- Sleep mode for extended battery life (Optional)
- Low battery indication.
- Overload protection.





SPECIFICATIONS

SPECIFICATION	GTA 200	GTA 300	GTA 600
Capacity	200g	300g	600g
Readability	0.01g	0.01g	0.01g
Linearity	+0.01	0.01	0.02
Pan Size	Liquid Crystal Display 13 mm	-	-
Display	Height	-	-
Weight	1.8Kg	-	-
Dimension	190 (W)x200 (D)x80 (H)mm	-	-
Power	9V/500mA AC/DC Adaptorm	-	-
Environment	0° C-40° C (32° F-104° F) RH <85%	-	-



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485
WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.acmasglobal.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57
Phase- IV Kundli, Sonapat, Haryana - 131028
Land Line No.: +91 - 11- 47619688,
Fax: +91-11-47619788
E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,
Sheung Wan, Hong Kong (PRC)
Tel.: 0086-13929598046 | 0086- 18922303099
E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,
Unit No 50, Straine- 1, Off- 502, Moscow, Russia
E-mail.: info@acmasindia.com
Tel.: 0049- 79592345 | Email: russia@acmasindia.com