


BALL **MILL**

MODEL NO. - ACM- 82303



BALL MILL

Model No. : 82303



BALL MILL

Model No. : 82303

Weiber **Ball Mills** are ideally used for mixing and grinding, electrically operated having capacity of few gms. to 1 Kg and 2 Kg. fitted with FHP motor. The unit drives with maximum speed of 80 rpm. The jar is made of aluminum and having steel balls of different sizes to be put in jar for mixing and grinding pesticides powder. Jar is interchangeable and easily removable, light in weight and compact unit working on 200 volts AC with cap. 1 Kg. (Jar made of stainless steel).

b) --- do ---- with cap.2 kg.(Jar made of S.S).



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485
WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.acmasglobal.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57
Phase- IV Kundli, Sonapat, Haryana - 131028
Land Line No.: +91 - 11- 47619688,
Fax: +91-11-47619788
E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,
Sheung Wan, Hong Kong (PRC)
Tel.: 0086-13929598046 | 0086- 18922303099
E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,
Unit No 50, Straine- 1, Off- 502, Moscow, Russia
E-mail.: info@acmasindia.com
Tel.: 0049- 79592345 | Email: russia@acmasindia.com