



Kahn (Reciprocating) Shaking Machine



Kahn (Reciprocating) Shaking Machine

It is general-purpose machine and found very useful for mixing or agitation ingredients, solution in flasks, bottles, beakers and Kahn racks etc. Upper tray covered with rubber floor and sides of the tray are slotted to take parallel and cross bars cover with rubber to hold flasks, bottles beakers, Kahn racks etc. The bars are tightened with washers and fly nuts can be adjusted as per the requirement. cycles single phase

The machine can accommodate Kahn racks, bottles or flasks of different sizes by adjusting the compartments with the help of rods supplied with the machine. Fitted with ¼ HP motor and having mechanical adjustments for adjusting the motion.

D.c. motor special system shaker to work on 220 volt a.c. for variable and smooth function speed control by speed regulator.

Kahn (Reciprocating) Shaking Machine





Number of oscillations fixed at 280 +/- 5 oscillations per minute. Supplied complete with top platform rods, ON/OFF switch, plug and cord.

Power Requirements :

To work on 200/220 volts, single phase, 50 Hz, AC supply.

Optional Accessories :

Timer 0-30 minutes.

PLATFORM SIZE :

- (a) 20" x 13" x 4"
- (c) 20" x 13" x 4"
- (d) 31" x 17" x 4"
- (e) 36" x 18" x 4"

Kahn (Reciprocating) Shaking Machine

Kahn (Reciprocating) Shaking Machine

OPTIONAL ACCESSORIES

- (a) DC MOTOR DRIVE with CONTINUOUSLY VARIABLE ELECTRONIC SPEED REGULATOR in lieu of mechanical speed control and A.C. motor

| Technical Matrix: | | |
|-------------------|--------|----------------------|
| PLATFORM SIZE | BODY | SPEED MOTOR CAPACITY |
| (a) 20"x13"x4" | Open | Fixed 1 /4 HP |
| (b) 20"x13"x4" | Closed | Variable 1 /4 HP |
| (C) 31"x17"x4" | Closed | Variable 1 /4 HP |
| (d) 36"x18"x4" | Closed | Variable 1 /2 HP |





We are having **Representatives in 55 Countries**

please visit our

International Dealers Section

to contact your nearest **ACMAS** Representative



TECHNOLOGY with **HUMAN TOUCH**



CE | An ISO 9001:2208 | ISO 14001:2008 | ISO 13485
WHO:GMP Products | GLP Compliant Products

www.acmasindia.com | www.scientificlaboratoryequipments.com | www.measuring-meters.com

ACMAS TECHNOCRACY PVT. LTD.

SALES OFFICE (INDIA)

312-313, Vardhman Capital Mall, L.S.C. 10,
GulabiBagh, Delhi-110052, INDIA
Tel: 0091-11-23646703, 23643054
(M) +91-9313971681 | Email: india@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/f Wing Tat Comm, Bldg 97, Bonham Strand East,
Sheung Wan, HONG KONG (PRC).
Tel: 0086-13929598046 | 0086-18922303099
Email: hk@acmasindia.com

SHOWROOM

141, Rai Industrial Estate, Rai, Sonapat,
Haryana-131029, INDIA
Tel: 0091-0-9312219738 (M) +91-9717741167
Email: info@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street, Ozerkovsky Embankment,
Unit No 50 Straine-1, Off-502, Moscow, RUSSIA
Tel: 0049-79592345 | Email: russia@acmasindia.com