



# TABLE TOP SIEVE SHAKER

Model No. : ACM-42306- U





# Table Top Sieve Shaker

Model No. : ACM-42306- U



# Table Top Sieve Shaker

Model No. : ACM-42306- U

An ideal equipment for normal day to day laboratory purpose, suitable for 8 12 Test Sieve (Soundless, Vibration type). Fitted with F.H.P Electric Motor. It can accommodate 6 Test sieves of 8 Dai with Lid and pan. Speed can be regulated through variable speed regulator to work on 220 volts A.C, speed regulator to work on 220 volts A.C.



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485

WHO: GMP Products | GLP Compliant Products

[www.acmasindia.com](http://www.acmasindia.com) | [www.cleanroom-equipments.com](http://www.cleanroom-equipments.com) | [www.test-chambers.com](http://www.test-chambers.com)

## **ACMAS** TECHNOLOGIES PVT. LTD.

### CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57  
Phase- IV Kundli, Sonapat, HR 131028  
Land Line No.: +91 - 11- 47619688,  
Fax: +91-11-47619788  
E-mail.: [info@acmasindia.com](mailto:info@acmasindia.com)

### SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,  
Sheung Wan, Hong Kong (PRC)  
Tel.: 0086-13929598046 | 0086- 18922303099  
E-mail.: [hk@acmasindia.com](mailto:hk@acmasindia.com)

### SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,  
Unit No 50, Straine- 1, Off- 502, Moscow, Russia  
E-mail.: [info@acmasindia.com](mailto:info@acmasindia.com)  
Tel.: 0049- 79592345 | Email: [russia@acmasindia.com](mailto:russia@acmasindia.com)