



WAGNER`S SHAKING MACHINE

MODEL NO.- ACM-42304-U





WAGNER`S SHAKING MACHINE

MODEL NO.- ACM-42304-U





WAGNER`S SHAKING MACHINE

MODEL NO.- ACM-42304-U



Fixed speed approx: 30 R.P.M with BELT SYSTEM. Fitted with H.P.
electric motor to work on 220 volts. A.C. (without bottles).

To accommodate 6 BOTTLES in two rows.

To accommodate 12 BOTTLES in two rows.



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485

WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.cleanroom-equipments.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57

Phase- IV Kundli, Sonapat, Haryana - 131028

Land Line No.: +91 - 11- 47619688,

Fax: +91-11-47619788

E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,

Sheung Wan, Hong Kong (PRC)

Tel.: 0086-13929598046 | 0086- 18922303099

E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,

Unit No 50, Straine- 1, Off- 502, Moscow, Russia

E-mail.: info@acmasindia.com

Tel.: 0049- 79592345 | Email: russia@acmasindia.com