



# WRIST ACTION SHAKING MACHINES

MODEL NO.- ACM-42302- U



# WRIST ACTION SHAKING MACHINES



D.c. motor special system shaker to work on 220 volt a.c. for variable and smooth function speed control by speed regulator.

It is general-purpose machine and found very useful for mixing or agitation ingredients, solution in flasks, bottles, beakers and Kahn racks etc. Upper tray covered with rubber floor and sides of the tray are slotted to take parallel and cross bars cover with rubber to hold flasks, bottles beakers, Kahn racks etc. The bars are tightened with washers and fly nuts can be adjusted as per the requirement. Cycles single phase.

The machine can accommodate Kahn racks, bottles or flasks of different sizes by adjusting the compartments with the help of rods supplied with the machine.



Fitted with  $\frac{1}{4}$  HP motor and having mechanical adjustments for adjusting the motion.

Number of oscillations fixed at 280 +/- 5 oscillations per minute. Supplied complete with top platform rods, ON/OFF switch, plug and cord.

### **Holding Capacity**

(a) To hold 4 flasks Cap. 100 ml. to 1000 ml. by side clamps

(b) To hold 8 flasks Cap. 100 ml. to 1000 ml. by side clamps

### **Optional Accessories**

Timer 0-30 minutes.



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485  
WHO: GMP Products | GLP Compliant Products

[www.acmasindia.com](http://www.acmasindia.com) | [www.acmasglobal.com](http://www.acmasglobal.com) | [www.test-chambers.com](http://www.test-chambers.com)

## **ACMAS** TECHNOLOGIES PVT. LTD.

### **CORPORATE OFFICE/ SHOWROOM (INDIA)**

Plot No. 352-353, Sector – 57  
Phase- IV Kundli, Sonapat, Haryana - 131028  
Land Line No.: +91 - 11- 47619688,  
Fax: +91-11-47619788  
E-mail.: [info@acmasindia.com](mailto:info@acmasindia.com)

### **SALES OFFICE (HONG KONG)**

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,  
Sheung Wan, Hong Kong (PRC)  
Tel.: 0086-13929598046 | 0086- 18922303099  
E-mail.: [hk@acmasindia.com](mailto:hk@acmasindia.com)

### **SALES OFFICE (RUSSIA)**

Inmed Trade Street Ozerkovsky Embankment,  
Unit No 50, Straine- 1, Off- 502, Moscow, Russia  
E-mail.: [info@acmasindia.com](mailto:info@acmasindia.com)  
Tel.: 0049- 79592345 | Email: [russia@acmasindia.com](mailto:russia@acmasindia.com)