



Wrist Action

Shaking Machine



Wrist Action Shaking Machine

Our Wrist-Action® Laboratory Shaker is a work horse in life science, biotechnology, pharmacology, chemistry, quality control, plant pathology, soil extractions, aphaeresis, toxicology, microbiology, bacteriology, food and pesticide analysis.

Features

- Rugged construction, with timer, rated for continuous use
- Shaking speed action is variable, simulating a snap wrist motion for a swirling effect, but the oscillations are constant at 380 to 390 per minute
- Degree of agitation, gentle to violent, is controlled by an adjustable lever with dial graduated from 1 to 10
- Power unit has constant-speed motor and built-in timer, range to 60 minutes in 1-minute intervals

Wrist Action Shaking Machine

Technical Details

Side sleeves have interchangeable side arm shafts 19 mm diameter and clamps. Erlenmeyer flasks from 125 to 1,000 mL capacity can be shaken at one time without need for balancing. Full torque of the 1/8 HP (93 W) motor is always applied to the oscillating mechanism regardless of speed setting. Side arm shafts take 2-place clamps with slip-on fittings which provide capacity of 4 to 16 flasks, suitable bottles, test tubes, etc., with minimum diameter of 10 mm. Platform for attachment to top takes 6 to 10 Erlenmeyer flasks up to 5,000 mL, also sputum bottles, etc. Platform clamps of rubber slide in parallel grooves. Power unit is enclosed in enameled sheet metal housing on rubber feet; overall dimensions, 343 x 389 x 298 mm high. It has a heavy cast iron base with rugged construction and two side arms for holding flasks. Fitted with a F.H.P. heavy duty motor with variable speed. Supplied with cord and plug to work on 220 / 230 volts A.C. supply.

Wrist Action Shaking Machine

CAPACITY

- (a) To hold 4 flasks
- (b) To hold 4 flasks and top platform to hold 4 flasks
- (c) To hold 8 flasks
- (d) To hold 8 flasks and top platform to hold 4 flasks

Wrist Action Shaking Machine



We are having **Representatives in 55 Countries**

please visit our

International Dealers Section

to contact your nearest **ACMAS** Representative



TECHNOLOGY with **HUMAN TOUCH**



An ISO 9001:2208 | ISO 14001:2008 | ISO 13485
WHO:GMP Products | GLP Compliant Products

www.acmasindia.com | www.scientificlaboratoryequipments.com | www.measuring-meters.com

ACMAS TECHNOCRACY PVT. LTD.

SALES OFFICE (INDIA)

312-313, Vardhman Capital Mall, L.S.C. 10,
GulabiBagh, Delhi-110052, INDIA
Tel: 0091-11-23646703, 23643054
(M) +91-9313971681 | Email: india@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/f Wing Tat Comm, Bldg 97, Bonham Strand East,
Sheung Wan, HONG KONG (PRC).
Tel: 0086-13929598046 | 0086-18922303099
Email: hk@acmasindia.com

SHOWROOM

141, Rai Industrial Estate, Rai, Sonapat,
Haryana-131029, INDIA
Tel: 0091-0-9312219738 (M) +91-9717741167
Email: info@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street, Ozerkovsky Embankment,
Unit No 50 Straine-1, Off-502, Moscow, RUSSIA
Tel: 0049-79592345 | Email: russia@acmasindia.com