




GLOSS METER

MODEL NO. - ACM- 52302-J



GLOSS METER

MODEL NO. - ACM- 52302-J



Weiber **Gloss Meter** is a compact, portable, sensitive & accurate instrument available in two models:- mains operated laboratory purpose and handy battery operated for the measurement of different geometry 20 deg, 60 deg & 85 deg.

The above instruments can be used for the measurement of gloss/brightness of paints, opaque glasses, plastics, ceramics, porcelains, paper, cloths, pigments, flours, powder etc.

This instrument is also useful to measurement of Brightness/whiteness of Plantation white sugar and measurement of reflectance of fabrics.

TECHNICAL SPECIFICATIONS

Range	0 – 99.9 gloss units
Angle heads	20, 60 & 85 degree
Accuracy	0.1 unit
Power Supply	9V battery, Pencils cells
Display	3-1/2 LED
Sensor	Silicon Photodiode
Light Source	LED Lamp
Standards	Black, White tile



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485
WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.acmasglobal.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57
Phase- IV Kundli, Sonapat, Haryana - 131028
Land Line No.: +91 - 11- 47619688,
Fax: +91-11-47619788
E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,
Sheung Wan, Hong Kong (PRC)
Tel.: 0086-13929598046 | 0086- 18922303099
E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,
Unit No 50, Straine- 1, Off- 502, Moscow, Russia
E-mail.: info@acmasindia.com
Tel.: 0049- 79592345 | Email: russia@acmasindia.com