



MODEL NO. - ACM- 22305-J

STEREOPIC ZOOM BINOCULAR **MICROSCOPE**

STEREOPSPIC ZOOM BINOCULAR **MICROSCOPE**



TECHNICAL SPECIFICATIONS

Observation head	Binocular head is inclined at 45° and can be rotated through 360° and can be looked in any position by using thumb screws. The interpupillary distance is adjustable between 55 to 75mm and right eye tube has dioptric focussing adjustment.
Objective	The microscope has a set of two standard objectives 2x & 4x. The low & high powered objectives are selected by rotating the nosepiece through 90°, having parfocal magnification.
Eyepieces	The standard extra wide field eyepieces WF 10X is provided
Magnification	Standard magnification 20x & 40x with flat field.
Image	Erect, three dimensional, crisp and sharp image.
Focusing	By rack and pinion mechanism with adjustable tension collar.
Illumination	With integral illuminator, two lamps are fitted to provide both transmitted & incident illumination. The lamps can be switched on separately or together. Complete with a glass plate having two stage clips.
Power Supply	22V-240V AC or 110V-120V AC, 50-60 Hz.
Accessories	Vinyl cover, dust cover, duster & cleaning brush.
Packing	In wooden storing cabinet with lock & key.
Optional Accessories	An array of different stands, stages and attachments designed to fit a variety of specific requirements. Choice of two objective combination 1x-2x or 1.5x-3x. spare eyepieces WF10X and dark field attachment



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485
WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.acmasglobal.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57
Phase- IV Kundli, Sonapat, Haryana - 131028
Land Line No.: +91 - 11- 47619688,
Fax: +91-11-47619788
E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,
Sheung Wan, Hong Kong (PRC)
Tel.: 0086-13929598046 | 0086- 18922303099
E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,
Unit No 50, Straine- 1, Off- 502, Moscow, Russia
E-mail.: info@acmasindia.com
Tel.: 0049- 79592345 | Email: russia@acmasindia.com