



# WEIBER INVERTED MICROSCOPE

(MODEL: ACM-22330-J)





# WEIBER INVERTED MICROSCOPE

(MODEL: ACM-22330-J)





## SPECIFICATIONS

<b>Viewing Head</b>	:	Trinocular Head Inclined at 30°
<b>Eyepiece</b>	:	Extra Wide Field Eyepiece EW 10X
<b>Objective</b>	:	LWD Infintive PLAN Objective 4X, LWD Infintive PLAN Objective 40X, Infintive PLAN PHASE Objective PH 10X, Infintive PLAN PHASE Objective PH 20X,
<b>Nosepiece</b>	:	Quintuple Nosepiece
<b>Condenser</b>	:	ELWD Condenser NA 0.3. LWD 72mm, (Without Condenser 150mm)
<b>Phase Annulus</b>	:	10X-20X Phase Annulus Plate
<b>Stage</b>	:	Plain Stage 160x250m, Glass Insert, Auxiliary Stages 70x180mm
<b>Focusing</b>	:	Coaxial Coarse and fine adjustment.
<b>Illumination</b>	:	Halogen Lamp 6V 30W
<b>Filter</b>	:	Blue, Green and Ground Glass
<b>Dimension</b>	:	48 x 41 x 67cm. Approx
<b>Weight</b>	:	12 Kg Approx.



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485  
WHO: GMP Products | GLP Compliant Products

[www.acmasindia.com](http://www.acmasindia.com) | [www.cleanroom-equipments.com](http://www.cleanroom-equipments.com) | [www.histopathology-lab-equipments.com](http://www.histopathology-lab-equipments.com)

## **ACMAS** TECHNOLOGIES PVT. LTD.

### **CORPORATE OFFICE/ SHOWROOM (INDIA)**

Plot No. 352-353, Sector – 57  
Phase- IV Kundli, Sonapat, Haryana - 131028  
Land Line No.: +91 - 11- 47619688,  
Fax: +91-11-47619788  
E-mail.: [info@acmasindia.com](mailto:info@acmasindia.com)

### **SALES OFFICE (HONG KONG)**

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,  
Sheung Wan, Hong Kong (PRC)  
Tel.: 0086-13929598046 | 0086- 18922303099  
E-mail.: [hk@acmasindia.com](mailto:hk@acmasindia.com)

### **SALES OFFICE (RUSSIA)**

Inmed Trade Street Ozerkovsky Embankment,  
Unit No 50, Straine- 1, Off- 502, Moscow, Russia  
E-mail.: [info@acmasindia.com](mailto:info@acmasindia.com)  
Tel.: 0049- 79592345 | Email: [russia@acmasindia.com](mailto:russia@acmasindia.com)