



TRINOCULAR MICROSCOPE

MODEL NO. - ACM- 22303-J



TRINOCULAR **MICROSCOPE**



TECHNICAL SPECIFICATIONS

Observation head	Trinocular Head is equipped with standard 45° inclined observation tubes, rotatable through 360° and can be fixed in any direction with a clamp screw. Interpupillary distance from 53 to 75 mm & diopter adjustment are provided. Right hand slide lever provides 100% light to binocular head or vertical photo tube.
Stand	One piece die cast arm and large size rectangular, heavy stable base with transformer/electrical fittings concealed inside.
Nosepiece	Precise quadruple revolving nosepiece with accurate centering and positive click stops for smooth accurate rotation.
Stage	built in graduated mechanical stage size 135x120 mm is controlled by convenient low co-axial positioned knobs (either side as per requirement) for easy and smooth scan of the specimen slide over entire range of 50 x 75 mm.
Focussing	By vertical movement of stage, coarse adjustment by rack pinion and fine movement by screw lever mechanism having a graduated knob
Illumination	6V-20W Halogen lamp fitted directly under the field lenses operated through solid state variable controlled transformer 220 or 110 volts. Mirror attachment while working in day light is also provided
Objectives	Achromatic 4x - 10x - 40x sL and 100 XSL oil immersion.
Eyepiece	Wide field WF10X and H5X or H6X or H15X (any two pairs)
Packing	In full mould Styrofoam / Thermocol container.
Accessories	Vinyl cover, duster, cleaning brush, filter & Operating Manual.
Optional Accessories	Dark-field condenser, phase contrast equipment, polarizing attachment, stage micrometer, micrometer discs, spare objectives, micro photographic equipment and CCTV, wooden storing cabinet etc



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485
WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.acmasglobal.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57
Phase- IV Kundli, Sonapat, Haryana - 131028
Land Line No.: +91 - 11- 47619688,
Fax: +91-11-47619788
E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,
Sheung Wan, Hong Kong (PRC)
Tel.: 0086-13929598046 | 0086- 18922303099
E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,
Unit No 50, Straine- 1, Off- 502, Moscow, Russia
E-mail.: info@acmasindia.com
Tel.: 0049- 79592345 | Email: russia@acmasindia.com