

Models:

ACM - 54006 - Q



# PARAFFIN BATH

Models:

- ABX-Series



# PARAFFIN BATH

The Paraffin Section Mounting Bath has been designed for ease of use and maximum safety. The top has generous space for resting the hands and a deep tapered bowl allowing slides to be used at any angle. The four rubber feet and the low level of the whole instrument provide stability. To

prevent the ingress of water spillage's the controls have been deeply recessed. A drop-on lid is provided to cover the bowl when not in use.

Models:

- ABX-Series



# PARAFFIN BATH

# PARAFFIN BATH

The bowl is made from stainless steel to offer maximum resistance against corrosion. Both the bowl and the upper casing have a matt black PTFE coating providing visual contrast against wax sections.

The lower casing has an easy-to-clean epoxy coating.

DIMENSIONS	Bowl: Diameter (top) 222mm x Depth 85mm External: Length 322mm x Depth 275mm x Height 120mm
TEMPERATURE RANGE	40° - 70°C (+/- 1.5°) at 30°C ambient
POWER SUPPLY/RATING	220 - 240V 50/60Hz AC 250 Watts
WEIGHT	3 Kg-3.5 Kg (Approx)



We are having **Representatives in 55 Countries**

please visit our

**International Dealers Section**

to contact your nearest **ACMAS** Representative



**TECHNOLOGY** with **HUMAN TOUCH**



An ISO 9001:2208 | ISO 14001:2008 | ISO 13485  
WHO:GMP Products | GLP Compliant Products

[www.acmasindia.com](http://www.acmasindia.com) | [www.scientificlaboratoryequipments.com](http://www.scientificlaboratoryequipments.com) | [www.measuring-meters.com](http://www.measuring-meters.com)

## **ACMAS** TECHNOCRACY PVT. LTD.

### **SALES OFFICE (INDIA)**

312-313, Yardhman Capital Mall, L.S.C. 10, GulabiBagh, Delhi-110052, INDIA

Tel: 0091-11-23646703, 23643054

(M) +91-9313971681 | Email: [india@acmasindia.com](mailto:india@acmasindia.com)

### **SALES OFFICE (HONG KONG)**

Unit D 28 11/f Wing Tat Comm, Bldg 97, Bonham Strand East, Sheung Wan, HONG KONG (PRC).

Tel: 0086-13929598046 | 0086-18922303099

Email: [hk@acmasindia.com](mailto:hk@acmasindia.com)

### **SHOWROOM**

141, Rai Industrial Estate, Rai, Sonapat, Haryana-131029, INDIA

Tel: 0091-0-9312219738 (M) +91-9717741167

Email: [info@acmasindia.com](mailto:info@acmasindia.com)

### **SALES OFFICE (RUSSIA)**

Inmed Trade Street, Ozerkovsky Embankment, Unit No 50 Straine-1, Off-502, Moscow, RUSSIA

Tel: 0049-79592345 | Email: [russia@acmasindia.com](mailto:russia@acmasindia.com)