



An ISO 9001:2208 | ISO 14001:2008 | ISO 13485
WHO:GMP Products | GLP Compliant Products

Economical Microprocessor Based **Dissolved Oxygen Meter**



(Acm-DO-820 Series)



Economical Microprocessor Based **Dissolved Oxygen Meter**

The Acm-MP-DDS - series economical microprocessor based dissolved oxygen meter are very accurate and economical DO meters. These microprocessor based DO meters come with single and two point calibration and data hold facility. The most unique feature of these DO meters are that they are provided with USB interface and can transfer the data to computer system in the laboratory and can be attached to various program applications for advanced studies.. The ACM-DO-820 series DO meters have automatic correction for barometric pressure and salinity. This equipment is unique for the use in college laboratories, research laboratories, water treatment and industrial application.

ACMAS TECHNOCRACY PVT. LTD.

www.measuring-meters.com | www.acmasindia.com | www.scientificlaboratoryequipments.com

Features:

- Microprocessor based
- Single or Two points automatic calibration
- Hold function
- Automatic temperature compensation
- Auto-Power Control
- Adjustable barometric pressure and salinity correction
- Automatic messages
- Expanded memory stores and recall up to 48 Data Sets
- USB interface and application software
- Large LCD display



Economical Microprocessor Based **Dissolved Oxygen Meter**



ACMAS TECHNOCRACY PVT. LTD.

www.measuring-meters.com | www.acmasindia.com | www.scientificlaboratoryequipments.com



Technical Features

DO Range	0 ~ 20mg/L(ppm)
Resolution	0.1mg
DO Accuracy	±0.5mg/L
% Saturation of Oxygen	0 ~ 200.0%
Accuracy	± 10.00%
Temperature Range	0.0 ~ 100.0°C
Temperature Accuracy	±0.1°C
Temperature Compensation	0 ~ 40°C
Salinity Correction	0 ~ 35g/L
Barometric Pressure Correction	450 ~ 850mmHg
Calibration Points	2 points
Memory	48 Data Sets
Output	USB
Power	1*9V Battery
Battery Life	100 hours
Dimensions	185(L)*mm88(W)mm*32(H)mm
Weight	300g



www.measuring-meters.com | www.acmasindia.com | www.scientificlaboratoryequipments.com

ACMAS TECHNOCRACY PVT. LTD.

SALES OFFICE (INDIA)

312-313, Vardhman Capital Mall, L.S.C. 10, GulabiBagh, Delhi-110052, INDIA

Tel: 0091-11-23646703, 23643054

(M) +91-9313971681 | Email: india@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/f Wing Tat Comm, Bldg 97, Bonham Strand East, Sheung Wan, HONG KONG (PRC).

Tel: 0086-13929598046 | 0086-18922303099

Email: hk@acmasindia.com

SHOWROOM

141, Rai Industrial Estate, Rai, Sonapat, Haryana-131029, INDIA

Tel: 0091-0-9312219738 (M) +91-9717741167

Email: info@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street, Ozerkovsky Embankment, Unit No 50 Straine-1, Off-502, Moscow, RUSSIA

Tel: 0049-79592345 | Email: russia@acmasindia.com