

ELECTRONIC DIGITAL CHEMICAL BALANCE



MODEL NO. - ACM- 78091- L

Electronic Digital **Chemical Balance**

FEATURES

- Large LCD display for easy view.
- Standard RS 232 interface.
- Parts Counting & percent weighing.
- Below balance weighty facility.
- Die cast aluminum design for long term stability & accurate result
- Various weighing units like ct, gm, oz, dwt, mom, GN.
- ISO-GLP compliance printouts (optional)



SPECIFICATIONS

SPECIFICATION	CTA 200	CTA 300	CTA 600
Capacity	200g	300g	600g
Readability	0.01g	0.01g	0.01g
Linearity	0.01	0.01	0.02
Pan Size	Liquid Crystal Display 13 mm		
Display	Height		
Weight	1.8Kg		
Dimension	190 (W)x200 (D)x80 (H)mm		
Power	9V/500mA AC/DC Adaptorm		
Environment	0° C-40° C (32° F-104° F) RH <85%		



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485
WHO: GMP Products | GLP Compliant Products

www.acmasindia.com | www.acmasglobal.com | www.test-chambers.com

ACMAS TECHNOLOGIES PVT. LTD.

CORPORATE OFFICE/ SHOWROOM (INDIA)

Plot No. 352-353, Sector – 57
Phase- IV Kundli, Sonapat, Haryana - 131028
Land Line No.: +91 - 11- 47619688,
Fax: +91-11-47619788
E-mail.: info@acmasindia.com

SALES OFFICE (HONG KONG)

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,
Sheung Wan, Hong Kong (PRC)
Tel.: 0086-13929598046 | 0086- 18922303099
E-mail.: hk@acmasindia.com

SALES OFFICE (RUSSIA)

Inmed Trade Street Ozerkovsky Embankment,
Unit No 50, Straine- 1, Off- 502, Moscow, Russia
E-mail.: info@acmasindia.com
Tel.: 0049- 79592345 | Email: russia@acmasindia.com